

Respitech Medical Inc.
250 Ranck Ave.
Lancaster, Pa. 17602
800-399-0250

Syringe Certification

Granted to:

Hospital
1 Hospital Avenue
Lancaster, Pa 12345

Serial Number of Certified Syringe: 123456789

Manufacturer of Syringe: **Hans Rudolph**

Actual Measured Volume in Liters: **3.000**

Test Conditions: Temperature: 24 C Barometric Pressure: 754.7 mmHg

June 20, 2013

Technician: **John Smith**

Signature

John Smith

The serial number of the reference syringe was 2010SY757

The stated reference volume was 2.9955 Liters

NIST traceable through New Mexico State Metrology Laboratory

Certificate # 5278